

Table S4 Primers used to amplify the NCU09969 locus in 500 bp intervals. These amplicons were used as the probe for the small RNA Northern blot. Sequences are listed in 5' to 3' direction.

	Forward	Reverse
1	ATGGGGAAGCAGGGAGAT	AAGCGATATGCCGACTG
2	GACGATGTGTTAGGCGCT	CTGGCGCCTTCTTTAAC
3	TACAGCAAAGAAGACAGAGAATCG	CTTTAGCATGGACAGACTCGTC
4	CGCATATTGGTCGTGTGC	TTGGTCCCATTCTGATTTCG
5	GCACGCGAAGGACGTCT	TCCATGATCTGTAGATGGTAACG
6	AATGCAACGATTTATGGCATC	AAATCTTAGTGCCCATCAATAAAG
7	CAAAACATGACCGTCGCC	TGATGATGTCGGGACATCG
8	GCATTTGTCGCATTCTCGA	CATGTGGTCGTCATGTAGTCG
9	GACCAAGCTGTTGCGACAT	GGTCGGTCGATGTGCTTAAT
10	TCAAACACCTCAATGTCCATTT	ATGTTGTGGTCGGATGTTCG
11	AAACTTCGACTGCAGTGTTC	GAAATTGATTGATTGAGACTAGTGG
12	TTTGTCTCTCGTCTCTCATT	GCTTCTGGTTGCGCTTCTAC