

Table S4 PCR Primers used to amplify Glyma15g08680 for CD-5.

PCR Reaction	Type	Primer Pair First Reaction	Type	Primer Pair Second Reaction (Nested)
1	Forward	CGGAGACTGGTAAATGTGAGC	Forward	CGGAGACTGGTAAATGTGAGC
1	Reverse	CAGCACACCTCCAAAACAAG	Reverse	CAGCACACCTCCAAAACAAG
2	Forward	GGCTAGGCCTTTGTGTTTGA	Forward	GGCTAGGCCTTTGTGTTTGA
2	Reverse	AACGGGAAATGCTGATTGAG	Reverse	AACGGGAAATGCTGATTGAG
3	Forward	GCTGCTGGATAGGTTTGAA	Forward	GCTGCTGGATAGGTTTGAA
3	Reverse	AACGGGAAATGCTGATTGAG	Reverse	AACGGGAAATGCTGATTGAG
4	Forward	TCAATTGCGGTTCTTCCTTC	Forward	ACGTCAGCTGCTATGAATGG
4	Reverse	TGGCCTCTGGCTTTCTTTAG	Reverse	GTAGCTCCAATGGCAAATC
5	Forward	TCAATTGCGGTTCTTCCTTC	Forward	TCAATTGCGGTTCTTCCTTC
5	Reverse	TGGCCTCTGGCTTTCTTTAG	Reverse	TGGCCTCTGGCTTTCTTTAG